

NEW

Remodeled Ophthalmoscope

BX α Plus

PURPOSE:

For evaluating and training Vision Fixation State and Preferred Retinal Locus (PRL)

For instructing patients on how to use their eyes in low-vision rehabilitation

ADVANTAGES:

Built-in various visual targets such as characters and figures

Using the absorptive filter (CCP-UG) therefore reducing glare 80%

A portable terminal such as iPhone[®] can be attached as a recording device (under development)




**Wow!
Amazing!
I can see!**



No Filter

CCP-UG Filter



The visual target "bird" is also available.